

Ruhr-Universität Bochum – Distribution of CLINITEST® Rapid-COVID-19-Antigen Self-Test:

This test product is approved by the BfArM as a lay test and is delivered by post in original packages containing 20 tests.


The original packages are to be handed out to two employees who will share the package. Wash hands thoroughly beforehand and wear a medical face mask.

Test ration per staff member:


- 10 swabs (1),
- 10 test cassettes (3),
- 10 test tubes with 10 pipette caps (4) (place in the plastic bag provided)
- 1 bottle with extraction buffer (2) (*an open bottle remains stable for up to 2 years if stored properly (2–30°C, closed cap)*).
- QR code for instructions



The cardboard workstation (5) should be split up.



Splitting the workstation



Disposal:

For individual tests: the used test kits are to be disposed of at home in a sturdy, tear-proof and sealed rubbish bag in the non-recyclable waste container. For (exceptional) group tests at RUB: the test kits are to be disposed of in sealed tear-proof and thick-walled rubbish bags in the non-recyclable waste bins outside the MA building.